

HMP_MDG_454DEFAULTPROTOCOL.V4.2

1. PURPOSE

- 1.1.** The current outline of the pilot production study is given below. This study is designed to meet a number of objectives including:
- 1.1.1. Test the use of the provisional 454 16S protocol thus providing additional feedback on the use of this protocol to generate 16S reads with actual samples
 - 1.1.2. Further inform quality control procedures of 16S reads generated from 454 beyond what can be learned with the Mock Community
 - 1.1.3. Further inform data analysis procedures of 16S reads generated from 454 beyond what can be learned with the Mock Community
 - 1.1.4. Allow centers to test in house LIMs and workflows and use this information to craft efficient production procedures
 - 1.1.5. Test procedures for data submission with the DACC and NCBI
 - 1.1.6. Test screening procedures for human contamination

2. IMPORTANT PROTOCOL NOTES

- 2.1.** This protocol describes the protocol for the clinical sample pilot study using barcoded primers for 16S variable regions V1-3 and V3-5.
- 2.1.1. DNA provided by Baylor for Pilot Experiment: clinical samples and HMP even mock minus Candida (positive control)
 - 2.1.1.1. Centers will quantify the samples using a fluorescent based assay and record this information.
 - 2.1.1.2. 1ul of each undiluted sample will be used per amplification reaction – samples may be diluted for accommodate robotic pipetting capacity, however the final concentration must be equivalent to the 1ul undiluted sample.
 - 2.1.1.3. Two controls will be amplified:
 - 2.1.1.3.1. The HMP even mock (1ul/amplification reaction) - positive control.
 - 2.1.1.3.2. One barcoded primer pair will be tested with water as a negative control.
 - 2.1.2. Each sequencing center will receive clinical samples from 12 individuals (12 individuals x 16.5 samples/individual = 198 samples).
 - 2.1.2.1. Each sample will be amplified with 2, 16S variable regions: (198 samples x 2 variable regions = 396 amplicons)
 - 2.1.2.2. Each sample will have a unique barcode (using either Broad or JCVI primers – see Appendix below).
 - 2.1.2.3. The 396 PCR products that will be cleaned with AmPure bead (SPRI), quantified (Quant-iT, or SybrGreen), normalized and pooled:
 - 2.1.2.3.1. 96 (up to 100) uniquely barcoded amplicons/pool totaling 4 pools (2 pools for V1-3 amplicons, and 2 pools for V3-5 amplicons).

- 2.1.2.3.2. Each 96 library pool will be sequenced on 1 PTP (~10,000 reads/amplicon) totaling 4 PTPs.
- 2.1.2.3.3. Alternative: 48 (up to 50) uniquely barcoded amplicons/pool totaling 8 pools (4 pools for V1-3 and 4 pools for V3-5 amplicons).
- 2.1.2.3.4. Each 48 library pool will be sequenced on ½ PTP (10,000reads/amplicon) totaling 8, ½ PTPs.
- 2.1.2.3.5. The HMP mock positive control should be included in each amplicon pool (96 or 48). For this study we are using V 1.1.
- 2.1.2.3.6. A negative control (water blank) should be included with each PCR core mix.
- 2.1.2.3.7. Note- multiplexed sample pools should consist of only one V region per pool (V1-V3 or V3-V5) and that pool should be run on a distinct half or full plate (in other words do not mix V regions in a pool or on the same plate ‘real estate.’)
- 2.1.2.3.8. If there are additional samples that cannot be accommodated in the pooling strategy desired by the center those samples (up to 8 samples) can be excluded (**Note- this is still under discussion by the centers**)
- 2.1.2.3.9. If the number of successfully amplified samples exceeds what a center can pool using one of the above strategies using four PTPs then the extra samples can be excluded from multiplexing.
- 2.1.2.3.10. **PCR failures- (Note- this is still under discussion by the centers)** if a PCR for a particular sample fails then the reaction should be repeated one time using the provisional protocol.
 - 2.1.2.3.10.1. If the reaction fails a second time, it should be repeated with increased template (double the volume) if the sample was a low concentration sample (e.g., <50ng/ul) or diluted to one half the concentration if the concentration is high (e.g., >50ng/ul).’Plan B’
 - 2.1.2.3.10.2. If the reaction fails again, the reaction can be considered a failure. However centers are free to try other adjustments to the PCR protocol at their discretion.
 - 2.1.2.3.10.3. All deviations from the default PCR conditions should be noted in the metadata file.

3. REQUIREMENTS

Materials/Equipment	Vendor	Catalog Number
AccuPrime™ Taq DNA Polymerase High Fidelity	Invitrogen	12346-086
Forward and Reverse Primers premixed	Operon	custom order
96 well thermocycler plate		
clear adhesive plate seals		
DNase/RNase free water	-	-
Thermo Cycler	-	-
Vortex	-	-
Pipettes	-	-
Aerosol resistant pipette tips		
MinElute PCR Purification Kit	Qiagen	28004
Ampure (SPRI) Beads (60mL kit)	Agencourt	A29152
1x low TE, pH 8.0		
Quant-IT ds DNA Assay, high sensitivity	Invitrogen/Molecular Probes	Q33120

Documentation

Quant-iT ds DNA Assay protocol (manufacturer's specifications)
SybrGreen Assay protocol (manufacturer's specifications)
poolingCalculator.xls
MinElute PCR Purification Kit Manual

4. METHOD

The PCR will be carried out using AccuPrime Taq High Fidelity. It is not necessary to setup this reaction on ice, however it is recommended.

- **PCR Primer Setup** Set up of 10uM primer plates (combining barcoded A primer with non-barcoded B primer- see appendix below for primer & tag sequences):
 - 1:10 dilution of the 100uM stocks:
 - For each variable region, set up a working primer plate with 90ul of 1x low TE
 - Add 5ul of each barcoded primer A from 100uM plate to corresponding well position in 10uM dilution plate.
 - Add 5ul of the 100uM of corresponding B adapter to each well of the 96-well plate (final concentration 10uM primer pair).
 - Mix by pipetting up and down.
 - Working concentration of 4uM (2uM each primer)
 - Dilute the 10uM primers 1:2.5 in 1x low TE (add 150ul of TE to each well of the 10uM primer stock plate and mix

- Primers can be stamped out into multiple single use primer plates and store at -20oC until ready to use.

4.1. PCR Setup - Mastermix

4.1.1. MasterMix contains the following amounts per sample:

13.85uL	RNAse/DNAse free water
2uL	10X AccuPrime PCR Buffer II
0.15uL	Accuprime Taq Hifi

16uL Total Volume of master mix

- 4.1.2. Multiply all volumes above by the amount of reactions needed plus 10%.
- 4.1.3. Combine reagents in a 2mL micro centrifuge tube and vortex to mix completely. If more than 100 reactions are needed a 15mL tube should be used.
- 4.1.4. Using an automated pipette transfer 16uL of master mix into individual wells in the 96 well reaction plate.
- 4.1.5. Cover plate and spin in a centrifuge at 2000rpm to collect sample at the bottom of the wells.

4.2. PCR Setup

- 4.2.1. For the initial reaction, transfer 2uL diluted (1:1 in water) DNA sample into the respective reaction wells. If the amplification fails repeat with 2uL of undiluted sample (See 2.1.1.2.)
- 4.2.2. Transfer 2uL of barcoded primers from primer plate to corresponding wells in 96 well PCR plate.
- 4.2.3. Securely seal with clear adhesive plate seal and vortex plate vigorously.
- 4.2.4. Spin briefly at 2000 rpm in a centrifuge.
- 4.2.5. Place in thermo cycler and cycle as follows:
- | | | |
|-------------|---------|-------------|
| 95°C | 2 min | } 30 cycles |
| 95°C | 20 sec | |
| 50 or 56°C* | 30 sec | |
| 72°C | 5min | |
| 4° | forever | |
- * 56°C for V3-1, 50°C for V5-3**
- 4.2.6. Clean PCR products using Agencourt AmPure Beads (use Agencourt protocol → 1.8x volume beads (36ul beads) – follow manufacturer’s specifications.
- 4.2.7. Elute beads with 25ul 1x low TE, pH 8.0 and transfer to new 96 well plate.

4.3. PCR Gel Analysis – (E-gel alternative using 1ul of PCR product - faster) - we will know from the Quantification step below if we have product so this step is actually optional.

- 4.3.1. In a new reaction plate add 1uL PCR product to 1uL 6X loading dye
- 4.3.2. Cover, vortex to mix, briefly centrifuge to collect sample at the bottom of the well.
- 4.3.3. Prepare a 1% agarose 1X TAE gel with EtBr.
- 4.3.4. Load samples and run approximately 1 hour at 100V.
- 4.3.5. Capture gel image on gel-doc and retain for analysis.

4.4. PCR Product Quantification

- 4.4.1. Quantify PCR product using SYBR-Green Quantification or Quant-IT ds DNA high sensitivity assay according to the manufacturer's specifications.

4.5. PCR Pooling

- 4.5.1. Using values from the SYBR Green or Quant-IT quantification, calculate pooling amounts using the poolingCalculator.xls or according to the following formula:

$$\text{Amount (uL) of each sample} = ((\text{vol}/2) * (\text{min})) / \text{sampleconc}$$

where:

Vol = total volume of each sample

Min = concentration in ng/ul of the sample with the lowest concentration

Sampleconc = concentration in ng/uL of target sample

- 4.5.2. Pool samples using a minimum transfer volume of 1uL. If less than 1uL is called for, a dilution must be made. If using the poolingCalculator.xls this will be accounted for.
- 4.5.3. Using a Qiagen minElute column, purify the pool according to the manufacturer's protocol.

(The Broad normalizes by converting all concentrations to molecules/ul. Determine which sample has the lowest concentration and then dilute all other samples to the same concentration. Pool equal volume of each (5-10ul) sample and then concentrate using a Qiagen MinElute column (elution with 30uL, 1x low TE, pH 8.0).

4.6. Sample Transfer for 454 Library Completion

- 4.6.1. Proceed directly to the qPCR library step.

Optional: Enter emPCR using ¼ the recommended primer concentration to avoid too many molecules amplified on bead. This results in high signal intensities during run, which leads to higher mixed reads and shorter read lengths.

5. POST SEQUENCING- Data Transfer to the DACC

- 5.1. After sequencing, the raw *.sff (e.g., no deconvoluting of samples or screening for human contamination) will be deposited at the DACC ftp site. Centers are of course free to work with this data in house in parallel to the DACC efforts.

- 5.2.** A sample spreadsheet containing metadata and library construction information as set up by the DACC will be completed by all centers and will accompany the submission of the *.sff files
- 5.3.** The DACC will perform a screening for the presence of human sequence and report these results to the centers

APPENDIX:


BROAD INSTITUTE PRIMER SEQUENCES INCLUDING TAGS

Purpose: In this approach, we will use 454 barcode sequences designed internally by the Broad (Pablo Alvarez and Will Brockman) between the A adapter and primer specific sequence (see picture below). Barcoded primer sets have been tested by the Broad Institute.


454 Barcoding (in 96 sample batches)

Clinical samples (96)

+ 96 barcoded primers



↓ PCR in 96 well plate




↓ Clean up with AmPure beads – at this point, several wells can be selected randomly and PCR products run on E-gel to confirm product size

↓ Quantify with Quan-IT ds DNA high sensitivity kit

↓ Normalize/pool

↓ qPCR quantify, 454 sequence



Broad Primer Sequences:

Primers were obtained from Operon reconstituted in TE for concentration of 100uM (100pmol/ul)

Variable regions V3 --> V1

"B" adapter for XLR + V1-3 27F

5' CCTATCCCCTGCTGTCCTGGCAGTCTCAGAGAGTTTGATCCTGGCTCAG

plate position	barcode name	barcode sequence	Barcode city name	primer name	A barcoded adapter for XLR system + barcode + V1-3 534R primer
A1	v2bBar8L	CACGC	Kabul	XLR_534R_v2bBar8L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCAGCATTACCGCGGCTGCTGG
A2	v2bBar23L	CGCAAC	Tirana	XLR_534R_v2bBar23L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCAGCATTACCGCGGCTGCTGG
A3	v2bBar174L	TGAAGC	Algiers	XLR_534R_v2bBar174L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTGAAGCATTACCGCGGCTGCTGG
A4	v2bBar602L	ACTTGC	Canberra	XLR_534R_v2bBar602L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACTTGCATTACCGCGGCTGCTGG
A5	v2bBar212L	TCACAC	Vienna	XLR_534R_v2bBar212L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTACACATTACCGCGGCTGCTGG
A6	v2bBar25L	CGTGAC	Baku	XLR_534R_v2bBar25L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCGTGACATTACCGCGGCTGCTGG
A7	v2bBar622L	ACGCGC	Nassau	XLR_534R_v2bBar622L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACGCGCATTACCGCGGCTGCTGG
A8	v2bBar72L	CCTCTC	Bridgetown	XLR_534R_v2bBar72L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCCTCTCATTACCGCGGCTGCTGG
A9	v2bBar600L	ACTCAC	Minsk	XLR_534R_v2bBar600L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACTCACATTACCGCGGCTGCTGG
A10	v2bBar559L	AGACAC	Brussels	XLR_534R_v2bBar559L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACTCAGACTTACCGCGGCTGCTGG
A11	v2bBar31L	CGACTC	Sarajevo	XLR_534R_v2bBar31L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGGACTCATTACCGCGGCTGCTGG
A12	v2bBar551L	AGCTTC	Rio	XLR_534R_v2bBar551L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAGCTTCATTACCGCGGCTGCTGG
B1	v2bBar1149L	AAGCGCG	Sofia	XLR_534R_v2bBar1149L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAAAGCGCATTACCGCGGCTGCTGG
B2	v2bBar15L	CAAGAAC	Ottawa	XLR_534R_v2bBar15L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCAAGAACATTACCGCGGCTGCTGG
B3	v2bBar556L	AGTTGGC	Bangui	XLR_534R_v2bBar556L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAGTTGGCATTACCGCGGCTGCTGG
B4	v2bBar144L	TATCAAC	Santiago	XLR_534R_v2bBar144L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGATCAACATTACCGCGGCTGCTGG
B5	v2bBar575L	AGGCGGC	Beijing	XLR_534R_v2bBar575L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAGCGGCATTACCGCGGCTGCTGG
B6	v2bBar48L	CGGTATC	Bogota	XLR_534R_v2bBar48L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCGGTATCATTACCGCGGCTGCTGG
B7	v2bBar166L	TGACGAC	Kinshasa	XLR_534R_v2bBar166L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTACGACATTACCGCGGCTGCTGG
B8	v2bBar613L	ACAAGGC	Brazzaville	XLR_534R_v2bBar613L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACAAGGCATTACCGCGGCTGCTGG
B9	v2bBar560L	AGACCTC	Zagreb	XLR_534R_v2bBar560L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAGACCTCATTACCGCGGCTGCTGG
B10	v2bBar741L	ATACCAC	Havana	XLR_534R_v2bBar741L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTACACATTACCGCGGCTGCTGG
B11	v2bBar228L	TGCGGCG	Nicosia	XLR_534R_v2bBar228L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTGCAGCATTACCGCGGCTGCTGG
B12	v2bBar807L	ATCTTAC	Prague	XLR_534R_v2bBar807L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGATCTTACATTACCGCGGCTGCTGG
C1	v2bBar1273L	AACCAGC	Copenhagen	XLR_534R_v2bBar1273L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTACAGCATTACCGCGGCTGCTGG
C2	v2bBar441L	TTCGAGC	Djibouti	XLR_534R_v2bBar441L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTTCGAGCATTACCGCGGCTGCTGG
C3	v2bBar1174L	AAGGTGC	Quito	XLR_534R_v2bBar1174L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAAGGTGCATTACCGCGGCTGCTGG
C4	v2bBar209L	TCTTGGC	Cairo	XLR_534R_v2bBar209L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTCTGGCATTACCGCGGCTGCTGG
C5	v2bBar153L	TAATCTC	Suva	XLR_534R_v2bBar153L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTAATCTCATTACCGCGGCTGCTGG
C6	v2bBar213L	TCACCTC	Helsinki	XLR_534R_v2bBar213L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTCACTCATTACCGCGGCTGCTGG
C7	v2bBar298L	TCCGCTC	Paris	XLR_534R_v2bBar298L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTCCGCTCATTACCGCGGCTGCTGG
C8	v2bBar146L	TATTGAC	Berlin	XLR_534R_v2bBar146L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTATTGACATTACCGCGGCTGCTGG
C9	v2bBar554L	AGTCGAC	Accra	XLR_534R_v2bBar554L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAGTCCGACATTACCGCGGCTGCTGG
C10	v2bBar646L	ACGGCTC	Athens	XLR_534R_v2bBar646L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAGCTCCATTACCGCGGCTGCTGG
C11	v2bBar158L	TGCGTTC	Guatemala	XLR_534R_v2bBar158L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTGCCTCATTACCGCGGCTGCTGG
C12	v2bBar207L	TCTCGAC	Conakry	XLR_534R_v2bBar207L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTCTCGACATTACCGCGGCTGCTGG
D1	v2bBar77L	CCAGGAC	Bissau	XLR_534R_v2bBar77L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCCAGGACATTACCGCGGCTGCTGG
D2	v2bBar601L	ACTCCTC	Budapest	XLR_534R_v2bBar601L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACTCCTCATTACCGCGGCTGCTGG
D3	v2bBar481L	TTCCTGC	Jakarta	XLR_534R_v2bBar481L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTTCCTGCATTACCGCGGCTGCTGG
D4	v2bBar419L	TTCATAC	Tehran	XLR_534R_v2bBar419L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTTCATACATTACCGCGGCTGCTGG
D5	v2bBar26L	CGTGCTC	Baghdad	XLR_534R_v2bBar26L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCTGCTCATTACCGCGGCTGCTGG
D6	v2bBar1172L	AAGGCAC	Dublin	XLR_534R_v2bBar1172L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAAGGCATTACCGCGGCTGCTGG
D7	v2bBar1210L	AACAACCTC	Jerusalem	XLR_534R_v2bBar1210L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAAACAATTACCGCGGCTGCTGG
D8	v2bBar606L	ACACGGAC	Rome	XLR_534R_v2bBar606L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACAGCATTACCGCGGCTGCTGG
D9	v2bBar159L	TGCCGAAC	Kingston	XLR_534R_v2bBar159L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTCCGAACATTACCGCGGCTGCTGG
D10	v2bBar147L	TATTGCTC	Tokyo	XLR_534R_v2bBar147L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTATTGCTCATTACCGCGGCTGCTGG
D11	v2bBar141L	TAGGAATC	Amman	XLR_534R_v2bBar141L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTAGGAATCATTACCGCGGCTGCTGG
D12	v2bBar119L	CCGGCCAC	Nairobi	XLR_534R_v2bBar119L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCCGGCCATTACCGCGGCTGCTGG
E1	v2bBar1379L	AATGGTAC	Tarawa	XLR_534R_v2bBar1379L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAAATGGTACATTACCGCGGCTGCTGG
E2	v2bBar208L	TCTCCGCT	Pyongyang	XLR_534R_v2bBar208L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTCTCCGCTCATTACCGCGGCTGCTGG
E3	v2bBar1267L	AACCTGGC	Seoul	XLR_534R_v2bBar1267L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAACCTGGCATTACCGCGGCTGCTGG
E4	v2bBar637L	ACGAAGTC	Bishkek	XLR_534R_v2bBar637L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACGAAGTCATTACCGCGGCTGCTGG
E5	v2bBar435L	TTCGTGGC	Riga	XLR_534R_v2bBar435L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTTCGTGGCATTACCGCGGCTGCTGG
E6	v2bBar1202L	AACACAAC	Beirut	XLR_534R_v2bBar1202L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAACAACATTACCGCGGCTGCTGG
E7	v2bBar413L	TTCCTGAC	Maseru	XLR_534R_v2bBar413L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTTCCTGACATTACCGCGGCTGCTGG
E8	v2bBar289L	TCCAAGTC	Monrovia	XLR_534R_v2bBar289L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTCCAAGTCATTACCGCGGCTGCTGG
E9	v2bBar433L	TTCGCGAC	Tripoli	XLR_534R_v2bBar433L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTTCGCGACATTACCGCGGCTGCTGG
E10	v2bBar121L	CCGGTCCG	Vaduz	XLR_534R_v2bBar121L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCCGGTCCGATTACCGCGGCTGCTGG
E11	v2bBar669L	ACCTGAAC	Vilnius	XLR_534R_v2bBar669L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACTGAAACATTACCGCGGCTGCTGG
E12	v2bBar1156L	AAGAGTTC	Luxembourg	XLR_534R_v2bBar1156L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAAGAGTTCATTACCGCGGCTGCTGG
F1	v2bBar370L	TTGACAAC	Bamako	XLR_534R_v2bBar370L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTTGACAACATTACCGCGGCTGCTGG
F2	v2bBar281L	TCCAAGAC	Valletta	XLR_534R_v2bBar281L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTCCAGAACATTACCGCGGCTGCTGG
F3	v2bBar49L	CGGTCTTC	Kishinev	XLR_534R_v2bBar49L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGCGGTCTTCATTACCGCGGCTGCTGG
F4	v2bBar1173L	AAGGCCTC	Monaco	XLR_534R_v2bBar1173L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAAGGCCTCATTACCGCGGCTGCTGG
F5	v2bBar599L	ACTAATTC	Rabat	XLR_534R_v2bBar599L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACTAATTCATTACCGCGGCTGCTGG
F6	v2bBar167L	TGACCGTC	Maputo	XLR_534R_v2bBar167L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTACCGTCATTACCGCGGCTGCTGG
F7	v2bBar161L	TGTCGGAC	Kathmandu	XLR_534R_v2bBar161L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTGTCCGACATTACCGCGGCTGCTGG
F8	v2bBar580L	AGGTTGTC	Amsterdam	XLR_534R_v2bBar580L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGAGTTGTCATTACCGCGGCTGCTGG
F9	v2bBar629L	ACGAGAAC	Wellington	XLR_534R_v2bBar629L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGACGAGAACATTACCGCGGCTGCTGG
F10	v2bBar184L	TGGTGAAC	Managua	XLR_534R_v2bBar184L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTGGTGAACATTACCGCGGCTGCTGG
F11	v2bBar233L	TCGTGTGC	Abuja	XLR_534R_v2bBar233L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTGTGTGCATTACCGCGGCTGCTGG
F12	v2bBar364L	TTGTGTTTC	Oslo	XLR_534R_v2bBar364L	CCATCTCATCCCCTGCGTGTCTCCGACTCAGTTGTGTTTCATTACCGCGGCTGCTGG

G1	v2bBar78L	CCACGGTC	Muscat	XLR_534R_v2bBar78L	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCAGCGCTATTACCCGGCTGCTGG
G2	v2bBar393L	TTGGAGGC	Islamabad	XLR_534R_v2bBar393L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTGGAGGCATTACCCGGCTGCTGG
G3	v2bBar350L	TTATCGGC	Asuncion	XLR_534R_v2bBar350L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTATCGGCATTACCCGGCTGCTGG
G4	v2bBar1164L	AAGAAGAC	Lima	XLR_534R_v2bBar1164L	CCATCTCATCCCTGCGTGTCTCCGACTCAGAAAGAACATTACCCGGCTGCTGG
G5	v2bBar1196L	AACTGTTC	Manila	XLR_534R_v2bBar1196L	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTTCATTACCCGGCTGCTGG
G6	v2bBar411L	TTCTCAAC	Warsaw	XLR_534R_v2bBar411L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTCAACATTACCCGGCTGCTGG
G7	v2bBar6L	CTTCCTTC	Lisbon	XLR_534R_v2bBar6L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTTCATTACCCGGCTGCTGG
G8	v2bBar1031L	ATTGCTAC	Doha	XLR_534R_v2bBar1031L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTGACATTACCCGGCTGCTGG
G9	v2bBar76L	CCTTCCGC	Moscow	XLR_534R_v2bBar76L	CCATCTCATCCCTGCGTGTCTCCGACTCAGCTTCCTCCGATTACCCGGCTGCTGG
G10	v2bBar555L	AGTCCGTC	Kigali	XLR_534R_v2bBar555L	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCCGTCATTACCCGGCTGCTGG
G11	v2bBar378L	TTGAACTC	Riyadh	XLR_534R_v2bBar378L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTGAACTCATTACCCGGCTGCTGG
G12	v2bBar1225L	AACGAGGC	Dakar	XLR_534R_v2bBar1225L	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACGAGGCATTACCCGGCTGCTGG
H1	v2bBar99L	CCGTTCCAC	Belgrade	XLR_534R_v2bBar99L	CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGTTCCACATTACCCGGCTGCTGG
H2	v2bBar236L	TCGAGAAC	Singapore	XLR_534R_v2bBar236L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGAGAACATTACCCGGCTGCTGG
H3	v2bBar731L	ACCGGAAGC	Bratislava	XLR_534R_v2bBar731L	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCGGAAGCATTACCCGGCTGCTGG
H4	v2bBar628L	ACGTTCCAC	Mogadishu	XLR_534R_v2bBar628L	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTTCACATTACCCGGCTGCTGG
H5	v2bBar1250L	AACGGAGTC	Pretoria	XLR_534R_v2bBar1250L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGAGTTCATTACCCGGCTGCTGG
H6	v2bBar438L	TTCGTTATC	Madrid	XLR_534R_v2bBar438L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTTATCATTACCCGGCTGCTGG
H7	v2bBar693L	ACCGTAATC	Colombo	XLR_534R_v2bBar693L	CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCTTAATCATTACCCGGCTGCTGG
H8	v2bBar672L	ACCTTGCTC	Khartoum	XLR_534R_v2bBar672L	CCATCTCATCCCTGCGTGTCTCCGACTCAGACCTTGCTCATTACCCGGCTGCTGG
H9	v2bBar355L	TTAAGATTC	Stockholm	XLR_534R_v2bBar355L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTAAGATTCATTACCCGGCTGCTGG
H10	v2bBar187L	TGGTTGGTC	Bern	XLR_534R_v2bBar187L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTGGTTGGTTCATTACCCGGCTGCTGG
H11	v2bBar162L	TGTCGGGTC	Damascus	XLR_534R_v2bBar162L	CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGGTCATTACCCGGCTGCTGG
H12	v2bBar1292L	AACCGTGC	Taipei	XLR_534R_v2bBar1292L	CCATCTCATCCCTGCGTGTCTCCGACTCAGAACGTCATTACCCGGCTGCTGG

Variable regions V5 --> V3

"B" adapter for XLR + V3-5 357F

5' CCTATCCCTGTGTGCTTGGCAGTCTCAGCCTACGGGAGGCAGCAG

plate position	barcode name	barcode sequence	Barcode city name	primer name	"A" barcoded adapter for XLR system + barcode + V3-5 926R primer
A1	v2bBar8L	CACGC	Kabul	XLR_926R_v2bBar8L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCAGCCCGTCAATTCMTTTRAGT
A2	v2bBar23L	CGCAAC	Tirana	XLR_926R_v2bBar23L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCGCAACCCGTCATTCMTTTRAGT
A3	v2bBar174L	TGAAGC	Algiers	XLR_926R_v2bBar174L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTGAAGCCCGTCAATTCMTTTRAGT
A4	v2bBar602L	ACTTGC	Canberra	XLR_926R_v2bBar602L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACTGTCGCCGTCATTCMTTTRAGT
A5	v2bBar212L	TCACAC	Vienna	XLR_926R_v2bBar212L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTCACACCCGTCATTCMTTTRAGT
A6	v2bBar25L	CGTAC	Baku	XLR_926R_v2bBar25L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTACCCGTCATTCMTTTRAGT
A7	v2bBar622L	ACGGGC	Nassau	XLR_926R_v2bBar622L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACGGCCGTCATTCMTTTRAGT
A8	v2bBar72L	CCTCTC	Bridgetown	XLR_926R_v2bBar72L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCCTCTCCCGTCATTCMTTTRAGT
A9	v2bBar600L	ACTCAC	Minsk	XLR_926R_v2bBar600L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCACCCTCAATTCMTTTRAGT
A10	v2bBar559L	AGACAC	Brussels	XLR_926R_v2bBar559L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACCCGTCATTCMTTTRAGT
A11	v2bBar31L	CGACTC	Sarajevo	XLR_926R_v2bBar31L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCGACTCCCGTCATTCMTTTRAGT
A12	v2bBar551L	AGCTTC	Rio	XLR_926R_v2bBar551L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAGCTCCCGTCATTCMTTTRAGT
B1	v2bBar1149L	AAGCCGC	Sofia	XLR_926R_v2bBar1149L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAAAGCCCGTCATTCMTTTRAGT
B2	v2bBar15L	CAAGAAC	Ottawa	XLR_926R_v2bBar15L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCAAGAACCCGTCATTCMTTTRAGT
B3	v2bBar556L	AGTTGGC	Bangui	XLR_926R_v2bBar556L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTGTGGCCCGTCATTCMTTTRAGT
B4	v2bBar144L	TATCAAC	Santiago	XLR_926R_v2bBar144L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGATATCAACCCGTCATTCMTTTRAGT
B5	v2bBar575L	AGGCGGC	Beijing	XLR_926R_v2bBar575L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAGGGCCGTCATTCMTTTRAGT
B6	v2bBar48L	CGGTATC	Bogota	XLR_926R_v2bBar48L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCGGTATCCCGTCATTCMTTTRAGT
B7	v2bBar166L	TGACGAC	Kinshasa	XLR_926R_v2bBar166L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTACGACCCGTCATTCMTTTRAGT
B8	v2bBar613L	ACAAGGC	Brazzaville	XLR_926R_v2bBar613L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACAAGGCCGTCATTCMTTTRAGT
B9	v2bBar560L	AGACCTC	Zagreb	XLR_926R_v2bBar560L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAGACCTCCCGTCATTCMTTTRAGT
B10	v2bBar741L	ATAACCAC	Havana	XLR_926R_v2bBar741L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGATACACCCGTCATTCMTTTRAGT
B11	v2bBar228L	TCGGGCG	Nicosia	XLR_926R_v2bBar228L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGGCCCGTCATTCMTTTRAGT
B12	v2bBar807L	ATCTTAC	Prague	XLR_926R_v2bBar807L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGATCTTACCCGTCATTCMTTTRAGT
C1	v2bBar1273L	AACCCAGC	Copenhagen	XLR_926R_v2bBar1273L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCCCGTCATTCMTTTRAGT
C2	v2bBar441L	TTCGAGC	Djibouti	XLR_926R_v2bBar441L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGAGCCCGTCATTCMTTTRAGT
C3	v2bBar1174L	AAGGTGC	Quito	XLR_926R_v2bBar1174L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGGTGCCGTCATTCMTTTRAGT
C4	v2bBar209L	TCTTGGC	Cairo	XLR_926R_v2bBar209L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTTGGCCCGTCATTCMTTTRAGT
C5	v2bBar153L	TAATCTC	Suva	XLR_926R_v2bBar153L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTAATCTCCCGTCATTCMTTTRAGT
C6	v2bBar213L	TCACCTC	Helsinki	XLR_926R_v2bBar213L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTCACTCCCGTCATTCMTTTRAGT
C7	v2bBar298L	TCCGCTC	Paris	XLR_926R_v2bBar298L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCCGTCATTCMTTTRAGT
C8	v2bBar146L	TATTGAC	Berlin	XLR_926R_v2bBar146L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGATTTGACCCGTCATTCMTTTRAGT
C9	v2bBar554L	AGTCGAC	Accra	XLR_926R_v2bBar554L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCCGACCCGTCATTCMTTTRAGT
C10	v2bBar646L	ACGGCTC	Athens	XLR_926R_v2bBar646L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACGGCTCCCGTCATTCMTTTRAGT
C11	v2bBar158L	TGCGTTC	Guatemala	XLR_926R_v2bBar158L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGCTCCCGTCATTCMTTTRAGT
C12	v2bBar207L	TCTCGAC	Conakry	XLR_926R_v2bBar207L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTGACCCGTCATTCMTTTRAGT
D1	v2bBar77L	CCAGGAC	Bissau	XLR_926R_v2bBar77L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCCAGGACCCGTCATTCMTTTRAGT
D2	v2bBar601L	ACTCCTC	Budapest	XLR_926R_v2bBar601L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACTCCCGTCATTCMTTTRAGT
D3	v2bBar481L	TTCTGTC	Jakarta	XLR_926R_v2bBar481L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTGCCGTCATTCMTTTRAGT
D4	v2bBar419L	TTCATAc	Tehran	XLR_926R_v2bBar419L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCATACCCGTCATTCMTTTRAGT
D5	v2bBar26L	CGTCGTC	Baghdad	XLR_926R_v2bBar26L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCGTCCCGTCATTCMTTTRAGT
D6	v2bBar1172L	AAGGCAC	Dublin	XLR_926R_v2bBar1172L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAAAGCCCGTCATTCMTTTRAGT
D7	v2bBar1210L	AACAACCTC	Jerusalem	XLR_926R_v2bBar1210L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAAACTCCCGTCATTCMTTTRAGT
D8	v2bBar606L	ACACGGAC	Rome	XLR_926R_v2bBar606L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACACCCCGTCATTCMTTTRAGT
D9	v2bBar159L	TGCCGAAc	Kingston	XLR_926R_v2bBar159L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGCAACCCGTCATTCMTTTRAGT
D10	v2bBar147L	TATTCGTC	Tokyo	XLR_926R_v2bBar147L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTATTCGTCGTCATTCMTTTRAGT
D11	v2bBar141L	TAGGAATC	Amman	XLR_926R_v2bBar141L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTAGGAATCCCGTCATTCMTTTRAGT
D12	v2bBar119L	CCGGCCAC	Nairobi	XLR_926R_v2bBar119L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGGCCACCCGTCATTCMTTTRAGT

E1	v2bBar1379L	AATGGTAC	Tarawa	XLR_926R_v2bBar1379L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAATGGTACCCGTC AATTCMTTTRAGT
E2	v2bBar208L	TCTCCGTC	Pyongyang	XLR_926R_v2bBar208L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTCTCCGTC CCGTC AATTCMTTTRAGT
E3	v2bBar1267L	AACCTGGC	Seoul	XLR_926R_v2bBar1267L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCTGGCCCGTCAATTCMTTTRAGT
E4	v2bBar637L	ACGAAGTC	Bishkek	XLR_926R_v2bBar637L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACGAGTCCCGTCAATTCMTTTRAGT
E5	v2bBar435L	TTCGTGGC	Riga	XLR_926R_v2bBar435L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTTGCGCCGTC AATTCMTTTRAGT
E6	v2bBar1202L	AACACAAC	Beirut	XLR_926R_v2bBar1202L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAACACAACCCCGTCAATTCMTTTRAGT
E7	v2bBar413L	TTCTTGAC	Maseru	XLR_926R_v2bBar413L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTGACCCGTC AATTCMTTTRAGT
E8	v2bBar289L	TCCAAGTC	Monrovia	XLR_926R_v2bBar289L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCAGTCCCGTCAATTCMTTTRAGT
E9	v2bBar433L	TTCGCGAC	Tripoli	XLR_926R_v2bBar433L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGCGACCCGTC AATTCMTTTRAGT
E10	v2bBar121L	CCGGTCCG	Vaduz	XLR_926R_v2bBar121L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCCGGTCCGCCGTC AATTCMTTTRAGT
E11	v2bBar669L	ACCTGAAC	Vilnius	XLR_926R_v2bBar669L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACCTGAACCCGTC AATTCMTTTRAGT
E12	v2bBar1156L	AAGAGTTC	Luxembourg	XLR_926R_v2bBar1156L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAGTCCCGTCAATTCMTTTRAGT
F1	v2bBar370L	TTGACAAC	Bamako	XLR_926R_v2bBar370L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGACAACCCGTC AATTCMTTTRAGT
F2	v2bBar281L	TCCAGAAC	Valletta	XLR_926R_v2bBar281L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCAGAACCCGTC AATTCMTTTRAGT
F3	v2bBar49L	CGGTCTTC	Kishinev	XLR_926R_v2bBar49L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCGGTCTCCCGTCAATTCMTTTRAGT
F4	v2bBar1173L	AAGCACTC	Monaco	XLR_926R_v2bBar1173L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGCCCTCCCGTCAATTCMTTTRAGT
F5	v2bBar599L	ACTAATTC	Rabat	XLR_926R_v2bBar599L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACTAATTCGCCGTC AATTCMTTTRAGT
F6	v2bBar167L	TGACCGTC	Maputo	XLR_926R_v2bBar167L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTGACCGTCCCGTCAATTCMTTTRAGT
F7	v2bBar161L	TGTCGGAC	Kathmandu	XLR_926R_v2bBar161L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTCCGACCCGTC AATTCMTTTRAGT
F8	v2bBar580L	AGGTGTGC	Amsterdam	XLR_926R_v2bBar580L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTTGTCCCGTCAATTCMTTTRAGT
F9	v2bBar629L	ACGAGAAC	Wellington	XLR_926R_v2bBar629L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACGAGAACC CGTCAATTCMTTTRAGT
F10	v2bBar184L	TGGTGAAC	Managua	XLR_926R_v2bBar184L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTGGTGAACCCGTC AATTCMTTTRAGT
F11	v2bBar233L	TCCGTTGC	Abuja	XLR_926R_v2bBar233L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTTGCCGTC AATTCMTTTRAGT
F12	v2bBar364L	TTGTGTTC	Oslo	XLR_926R_v2bBar364L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGTGTCCCGTCAATTCMTTTRAGT
G1	v2bBar78L	CCACGGTC	Muscat	XLR_926R_v2bBar78L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCCACGGTCCCGTCAATTCMTTTRAGT
G2	v2bBar393L	TTGGAGGC	Islamabad	XLR_926R_v2bBar393L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGGAGGCCCGTCAATTCMTTTRAGT
G3	v2bBar350L	TTATCGGC	Asuncion	XLR_926R_v2bBar350L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTATCGGCCCGTCAATTCMTTTRAGT
G4	v2bBar1164L	AAGAAGAC	Lima	XLR_926R_v2bBar1164L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGAAGACCCGTC AATTCMTTTRAGT
G5	v2bBar1196L	AACTGTTC	Manila	XLR_926R_v2bBar1196L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAAGTTCGCCGTC AATTCMTTTRAGT
G6	v2bBar411L	TTCTCAAC	Warsaw	XLR_926R_v2bBar411L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCTCAACCCGTC AATTCMTTTRAGT
G7	v2bBar6L	CTTCCTTC	Lisbon	XLR_926R_v2bBar6L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCCTTCCCGTCAATTCMTTTRAGT
G8	v2bBar1031L	ATTCTGAC	Doha	XLR_926R_v2bBar1031L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGATTCGTACCCGTC AATTCMTTTRAGT
G9	v2bBar76L	CCTTCCGC	Moscow	XLR_926R_v2bBar76L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGCCCTCCGCCGTC AATTCMTTTRAGT
G10	v2bBar555L	AGTCCGTC	Kigali	XLR_926R_v2bBar555L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAGTCCGTC CCGTCAATTCMTTTRAGT
G11	v2bBar378L	TTGAAGTC	Riyadh	XLR_926R_v2bBar378L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTGAAGTCCCGTCAATTCMTTTRAGT
G12	v2bBar1225L	AACGAGGC	Dakar	XLR_926R_v2bBar1225L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCAGCCCGTCAATTCMTTTRAGT
H1	v2bBar99L	CCGTTCCAC	Belgrade	XLR_926R_v2bBar99L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACGTTCCACCCGTC AATTCMTTTRAGT
H2	v2bBar236L	TCCGAGAAC	Singapore	XLR_926R_v2bBar236L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTCCGAGAACCCGTC AATTCMTTTRAGT
H3	v2bBar731L	ACCGGAAGC	Bratislava	XLR_926R_v2bBar731L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACCGGAAGCCCGTCAATTCMTTTRAGT
H4	v2bBar628L	ACGTCCAC	Mogadishu	XLR_926R_v2bBar628L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACGTTCCACCCGTC AATTCMTTTRAGT
H5	v2bBar1250L	AACGGAGTC	Pretoria	XLR_926R_v2bBar1250L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCAGGAGTCCCGTCAATTCMTTTRAGT
H6	v2bBar438L	TTCGTTATC	Madrid	XLR_926R_v2bBar438L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTCGTTATCCCGTCAATTCMTTTRAGT
H7	v2bBar693L	ACCGTAATC	Colombo	XLR_926R_v2bBar693L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACCGTAATTCGCCGTC AATTCMTTTRAGT
H8	v2bBar672L	ACCTTGGTTC	Khartoum	XLR_926R_v2bBar672L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGACCTTGGTCCCGTCAATTCMTTTRAGT
H9	v2bBar355L	TTAAGATTC	Stockholm	XLR_926R_v2bBar355L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTTAAGATTC CCGTCAATTCMTTTRAGT
H10	v2bBar187L	TGGTTGGTC	Bern	XLR_926R_v2bBar187L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTGGTTGGTCCCGTCAATTCMTTTRAGT
H11	v2bBar162L	TGTCGGGTC	Damascus	XLR_926R_v2bBar162L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGTGTCCGGTCCCGTCAATTCMTTTRAGT
H12	v2bBar1292L	AACCGTGTTC	Taipei	XLR_926R_v2bBar1292L	5' CCATCTCATCCCTGCGTGTCTCCGACTCAGAACCAGTCCCGTCAATTCMTTTRAGT

JCVI primers:

Please note that the 27F (V1-> V3 region) and the 357F (V3->V5) primers consist of 'B' Adapter and the respective 16S primer. Therefore they are the same sequence for each respective V region (you only need to order one 'F' primer to pair with each individual 'R' primer). Each 'R' primer (534 R (V1->V3) and 926R (V3->V5) listed contains a unique barcodes. The barcodes have also been posted in a tab-delimited file on the DACC site (jcvl_barcode_v1v3_v3v5_101409). Thanks to Kelvin Li for work on the barcode design.

